



At-line/on-line measurement of sugars in HMO fermentation

Katrine Lund Andersen, Krist V. Gernaey, Thomas Skov, Kathrine Bisgaard Christensen



Introduction

Health benefits of Human Milk Oligosaccharides



Immune system

Microbiome



Cognition and brain development

Bacterial balance



Intestinal barrier

Project goal: To evaluate the potential of process analytical technology to generate real-time quantification of key substrates and products in HMO fermentation as a first step towards development of future process control applications

Methods



2-FUCOSYLLACTOSE (2-FL)



3-FUCOSYLLACTOSE (3-FL)



3-SIALYLLACTOSE (3-SL)



LACTO-N-TETRAOSE (LNT)

Sample preparation


Vibrational spectroscopy

Multivariate data analysis


Model validation

Outcome

Successful measurements of components of interest



Models with R² and RMSEcv values indicating satisfactory accuracy of HMO quantification



Perspectives

- Investigation of in-line performance of the spectroscopic method
- Assessment of the robustness of the models to ensure optimal production line performance

A Helix Lab Fellowship

- Unique industry insights
- Great networking opportunity
- Resourceful support system and friends for life

